



Mill specific LIMS solution

Safe Computing delivers a Mill specific Laboratory Information Management Solution (LIMS).

Due to the integration between the LIMS module and the weighbridge it is possible to ensure that ingredients are held at the weighbridge awaiting positive release from QC upon arrival on-site at the mill.

Tests supported and recorded in the Safe LIMS system via integration to the equipment in the lab include the capture of wheat and flour analysis, such as moisture levels, wheat variety, specific gravity and screening percentages.

The Safe system will also provide support for the calculation of any claim due to be charged for non-conformance to contract, and will pass this information seamlessly back to the weighbridge and accounts modules for further processing if the delivery falls outside pre-defined tolerances against specification.

Production batches are also held awaiting QC if required and allow for samples to be tested before the system allows release to distribution.

Outbound conformance certificates of analysis are produced from the system to go with the goods whether bulk or bags or to be distributed electronically to the customer from the system.

